



PHILIPS

Customer Services
Clinical Education

Functional CT - Liver Perfusion

Purpose

Functional CT is intended for visualization, assessment and quantification of blood flow, blood volume, time to peak and peak enhancement using dynamic CT data. The application provides whole-organ or single-location liver, lung and kidney perfusion calculations.

Before you begin





Intravenous contrast bolus is injected into the patient and region of interest is scanned repeatedly for a period-of-time. Position the slice location to include the liver, spleen, portal vein, and aorta.

Inject Contrast. (Recommended injection parameters: 40-50 ml total, >300 mgI/ml, 5-7 ml/ second.)

Obtain axial, contrast-enhanced scans at the same location during patient's breath-hold for the duration, from the injection to the end of the portal phase. (Recommended scan parameters: 5 mm slice thickness, 1 second scan duration, 3 second cycle time, 360 degree rotation, 120kV, 150-200 mAs, breath hold.)

Workflow

To analyze a CT liver perfusion scan:


1. Select the study and series and click .
2. Click  (**Draw Artery**) and draw an arterial ellipse.
3. Click  (**Draw Vein**) and draw a venous ellipse.
4. In the **Perfusion Type** field, click the drop-down arrow and select **Liver**.
5. Click  (**Draw Spleen**) and draw a spleen ellipse.



The arterial, venous, and spleen ellipses do not need to be on the same image.

Each of the drawn ellipses are now labeled within the main viewport.








6. Click  next to Vessel Definition to move to step 2, **Perfusion Maps**. The perfusion maps display in the top right viewport.



Liver Perfusion maps include the following:

- Arterial Perfusion
- Portal Venous Perfusion
- Total Perfusion
- Hepatic Perfusion Index

7. Under **Tissue ROIs**, click  (**Draw tissue ROI/Ellipse**) as needed.
8. Do one or more of the following:
 - To save the images, click  (**Save**) and select the way you want to save the images.
 - To bookmark the images, click  (**Bookmark**) and set the bookmark in the **Save Bookmark** dialog box.
 - To send the final images to the reporting package, click the drop-down arrow next to the displayed reporting option  and select .

© 2020 Koninklijke Philips N.V. All Rights Reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without prior written consent of the copyright owner.

Limitations of Liability

Philips has taken care to ensure the accuracy of these materials. However, Philips assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Philips may make improvements or changes in the product(s) or program(s) described in this content at any time.

Educational use only

Philips is providing these materials for educational use only. The content is not intended to provide medical advice, diagnosis, or treatment. This content should not be regarded as a substitute or replacement for the Operating Manual or Instructions for Use. Before attempting to use any Philips medical devices, you must read the Instructions for Use thoroughly, paying particular attention to all warnings, cautions, and notes it may contain. In all cases the information contained in the Instructions for Use is considered leading information.

Terms of use

Please see the full terms of use. [HTTP://philips.com/terms](http://philips.com/terms)

