



NeuroQ Create Database



Purpose

This quick step provides an optimal workflow scenario. Please keep in mind this application has additional workflow steps and functionality, which are detailed in the *Instructions for Use*.



Before you begin

- Format a Flash Drive for ISP.
- Create a minimum of two NeuroQ Results Files.

Workflow

To create a NeuroQ Database:

1. Insert a flash drive in the computer.
2. Launch NeuroQ and process a study from an asymptomatic control patient. Select **Save** to create a NeuroQ Review file. Run NeuroQ additional times to create review files from multiple control patients.

NOTE: You must have at least two review files available in order to create a new database.

3. Copy the NeuroQ review files (which have a description of NeuroQ Review File Reformatted:xx) from the study to the flash drive using the IntelliSpace Portal Patient Save Directory Copy function. The review files are DICOM secondary capture files saved in the Study\Series directory.
4. Select a valid series and launch NeuroQ. The study is not needed for the database creation, but is only needed for launching NeuroQ successfully.
5. Select **Control Database > Create New Database**. The Select Files window opens.
6. Select the review files that you need to include in a new database.
7. Type a name for the database file.
8. Select **Create DB**.
9. Once you create the database, it is available for use in subsequent NeuroQ sessions when you select **Control Database > Select Database** within NeuroQ and pick the folder name matching the database created above.

© 2021 Koninklijke Philips N.V. All Rights Reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without prior written consent of the copyright owner.

Limitations of Liability

Philips has taken care to ensure the accuracy of these materials. However, Philips assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Philips may make improvements or changes in the product(s) or program(s) described in this content at any time.

Educational use only

Philips is providing these materials for educational use only. The content is not intended to provide medical advice, diagnosis, or treatment. This content should not be regarded as a substitute or replacement for the Operating Manual or Instructions for Use. Before attempting to use any Philips medical devices, you must read the Instructions for Use thoroughly, paying particular attention to all warnings, cautions, and notes it may contain. In all cases the information contained in the Instructions for Use is considered leading information.

Terms of use

Please see the full terms of use. [HTTP://philips.com/terms](http://philips.com/terms)

